



Fastback Protein A Sepharose® FF Resin Datasheet

Fastback Protein A Sepharose® FF Resin designed for simple, one-step and rapid antibody purification from serum, ascites and tissue culture supernatant such as those derived from static cultures and bioreactors. Recombinant protein A has been coupled to Sepharose® to obtain a stable matrix with a high binding capacity for immunoglobulins via the heavy chain of the FC region (up to 30 mg/ml Human IgG). Antibody samples purified using this affinity resin may be used in a wide range of laboratory procedures such as 1D or 2D polyacrylamide gel electrophoresis, Western blotting, ELISA etc. Binding affinity varies depending upon the source species and subclass see our website for further details.

Specification:

Specificity:	Protein A affinity antibodies
Matrix:	Sepharose®
Coupled ligand:	3.5 mg Protein A/ml resin
Binding capacity (human IgG):	30 mg/ml
Bead size:	60-165 µm
Flow rate:	0.25-1 ml/min (recommended)
Maximum pressure:	120-140 psi
Buffer compatibility:	Common aqueous buffers from pH 2.5-10
Binding buffer example:	1.5 M glycine/NaOH, 3 M NaCl pH 9.0
Elution buffer example (1):	0.1 M sodium citrate pH 5.5
Elution buffer example (2):	0.2 M glycine/HCl pH 2.5
Neutralization buffer example:	1 M tris/HCl pH 9.0
Shipping/delivery:	50% (v/v) resin suspension in 0.01% thimerosal
Storage:	0.01% thimerosal at 2-8°C for up to 2 years from manufacture

Ordering Information:

Product	Volume	Order Code
Fastback Protein A Sepharose® FF Resin (1 ml)	1 ml	Fastback-PA-1
Fastback Protein A Sepharose® FF Resin (5 ml)	5 ml	Fastback-PA-5
Fastback Protein A Sepharose® FF Resin (25 ml)	25 ml	Fastback-PA-25

Protein Ark Limited

Telephone +44 (0) 114 224 2236
FAX: +44 (0) 114 224 2222
Email: info@proteinark.com



The Innovation Centre, 217 Portobello, Sheffield S1 4DP, UK.
Email: info@proteinark.com Telephone: +44 (0) 114 224 2236
Web: www.proteinark.com Fax: +44 (0) 114 224 2222