

## Proteus Mini Purification Spin Column Pack

Protein Ark's Proteus mini purification spin columns are designed for small scale protein purifications such as those required for expression trials, solubility determination tests, screening, titering and scouting studies. The columns fit all standard microfuges and allow you to process multiple samples in parallel.



### Specification:

Sinter type:	Ultra High Density Polyethylene
Plastic construction:	Polypropylene
Pore size:	10-40 $\mu\text{m}$ (mean pore size, 18 $\mu\text{m}$ )
Maximum sample volume:	0.65 ml
Hold-up volume	< 5 $\mu\text{l}$
Maximum g force:	14,000 g
Typical spin times:	1 min for 0.65 ml sample at 14,000 g
Storage:	Store at RT (non-sterile)
Shelf-life:	24 months

### Ordering Information:

Product	Units	Order Code
Proteus Mini Purification Spin Column Pack	25	GEN-MPS500

### Protein Ark Limited

Telephone +44 (0) 33 33 44 20 25  
FAX: +44 (0) 33 33 44 20 25  
Email: [info@proteinark.com](mailto:info@proteinark.com)