

Realtime stain

The only Universal Protein Pre-stain on the Market

Make your own in-house protein ladders (Lane 1)

View protein bands as they run (Lanes 3-12)

Behaves like Coomassie


Quick 'look & see'

'WOW'

'Amazing'

'Cool'

New Kid on the Block



Visualise your samples as they run on your SDS gel

Example Gels: NuSep Tris-Glycine 4-20% SDS PAGE. Run using standard method and NEW Realtime stain. Images taken in cassette directly after running.

| Description: | No. of gels (10 lanes) | No. of lanes | Order codes: |
|------------------------|------------------------|--------------|-------------------|
| Realtime Stain (2ml) | 40 | 400 | GEN-RT-Stain-2000 |
| Realtime Stain (0.2ml) | 4 | 40 | GEN-RT-Stain-200 |