

Fact Sheet

Eco Pre-stained Protein Ladder

3 COLOR PRE-STAINED PROTEIN LADDER

HIGHLY CONCENTRATED

ECONOMICAL

Features

- 9 Pre-Stained Proteins
- 3 Colors - No more guessing
- Supplied in gel-loading buffer, ready to use
- Wide range molecular weight 15-180 kDa
- Store at -20°C, up to 24 months

Designed for


- SDS-polyacrylamide gels
- Western blots (any membrane: PVDF, nylon, nitrocellulose)

Cost Analysis

- Highly Concentrated - less usage per lane
- Up to 250 Gels per 2x250µL package

Ordering Information

| Product Name | Product Code | Price |
|---------------------------------------|--------------|--------|
| Eco Pre-stained Protein Ladder, 500µL | PAL-ECL-500 | 80 GBP |



Use with Quick Coomassie Gel Stain for even more savings!

Price per gel stained: **£1.00**
(20 mL/gel, used 3x)

Time to stain a gel: **15 min**

Save money and time!

